

FoLix™ 



TRANSFORMING
HAIR LOSS
TREATMENTS
WITH LASER

LUMENIS

LUMENIS.COM

CAPTURE THE DEMAND




2 in 3 men
face hair loss
by the age of 35¹

1 in 2 women
face hair loss
by the age of 50²

66% of women think about their hair loss once a day or more.³

44% of men between 16 and 24 years are 'worried' about hair loss.⁴



¹ American Hair Loss Association

² MedlinePlus by NIH

³ Viviscal survey - Rahalhairtransplant.com

⁴ Sanghvi, S. (2022) UK male hair loss report. LloydsPharmacy

⁵ As measured in 1-month and 3-months follow ups following last treatment session, compared to baseline. Data on file.

GROW WITH FoLix

FoLix is a revolutionary non-ablative fractional laser system using hair stimulation technology, providing effective, safe, and simple treatment for women and men.



First FDA-cleared

The first FDA-cleared fractional laser for hair loss treatment with FLX™ technology.



Natural and safe - No chemicals

Leveraging the body's own processes for hair stimulation without chemicals.



Proven results

Clinical studies show increase in hair count and improvement in scalp hair appearance.⁵



Simple, 30 min procedure

Simple to use, easy to learn, quick to treat.



No anesthesia. No needles. No downtime.




TRUST SCIENCE TRUST FoLix



Unlock Nature with Science

FoLix uses targeted laser energy for dermal coagulation, activating natural healing processes. This boosts blood flow, cytokine activity, and hair follicle stimulation.



Histology shows the columns of coagulative zones created after laser*

Clinical Studies Show**

10-18% increase in number of hairs per cm².

13-16% additional follicular units per cm².



"People are looking for non-surgical, non-pharmaceutical solutions, and that's where FoLix comes in."

Marc Avram, MD

Harvard trained, board-certified dermatologist who specializes in hair loss and hair transplantation. Clinical professor of dermatology at Cornell Medical Center.







"FoLix represents an opportunity to differentiate your practice with effective laser innovation from a global leader in medical aesthetics."

Neil Sadick, MD

A leading recognized dermatologist involved in the research of advanced technologies and treatments for hair loss. Clinical Professor of Dermatology at Weill Cornell Medical College.

CLINICALLY PROVEN RESULTS



Courtesy of LaserMed clinic

*Taken on ex-vivo porcine tissue

**As measured in 1-month and 3-months follow ups following last treatment session, compared to baseline. Data on file.

SIMPLE TO USE EASY TO LEARN QUICK TO TREAT



Fast, 30 minutes treatment with no anesthesia and minimal discomfort.



3-6 treatment sessions to achieve results.



4 weeks intervals plus maintenance 2-3 times a year.

Intuitive User Interface

Incorporates FDA-cleared settings to ensure efficient, consistent treatments.

FoliCool Tip (Patent Pending)

Enhances safety and comfort with constant cooling. Designed to separate the hair for smooth treatment.

FLX™ Laser Technology

Unique algorithm ensures optimal irradiance and area coverage to maximize operational efficiency, patient comfort, and clinical outcomes.



About Lumenis:

Lumenis is a global leader in the medical aesthetic market and is a world-renowned expert in developing and commercializing innovative energy-based technologies, including Laser, Intense Pulsed Light (IPL) and Radio-Frequency (RF). For over 50 years, Lumenis' ground-breaking products have redefined medical treatments and set technological and clinical gold-standards, revolutionizing existing treatment methods and creating solutions for previously untreatable conditions.

Efficient system support - advanced remote connectivity for swift diagnostics and resolutions of service issues, maintaining clinic operations without disruption.

Indication for use:

Improving the appearance of scalp hair in adult males and females with Fitzpatrick skin types I to IV, who are seeking treatment for hair loss.

Warnings and risks:

FoLix treatment could cause redness, swelling, scarring, and pigmentation change. The use of FoLix is contraindicated for patients with any concurrent cancer or history of skin cancer, active infection, fungal or bacterial diseases. See the system user manual for a complete list of contraindications and risks.

Testimonial represents a specific experience of a paid spokesperson. Experience of others may vary.



© 2024 All Rights Reserved.
The Lumenis Be Group of Companies.

Manufactured for Lumenis Be Ltd.
Yokneam Industrial Park
Hakidma Street 9, P.O.B. # 426
Yokneam 2069236, Israel
Tel: +972.77.959.9000
www.Lumenis.com

AMERICAS
San Jose, CA, USA
T +1 408 764 3000
T +1 877 586 3647
F +1 408 764 3999

PB-00076330 Rev A